

```
1 // Annex II
2 // Main window for conducting and analysing the tremor experiments with the g.Nautilus accelerometer
3 // Written by Javi Rodriguez and Ben Hesketh
4
5 using System;
6 using System.Collections.Generic;
7 using System.Linq;
8 using System.Text;
9 using System.Threading.Tasks;
10 using System.Windows;
11 using System.Windows.Controls;
12 using System.Windows.Data;
13 using System.Windows.Documents;
14 using System.Windows.Input;
15 using System.Windows.Media;
16 using System.Windows.Media.Imaging;
17 using System.Windows.Navigation;
18 using System.Windows.Shapes;
19 using Gtec.Gds.Client.API.Wrapper;
20 using Gtec2;
21 using System.Timers;
22 using TFM;
23
24 namespace accelerometer
25 {
26     /// <summary>
27     /// Interaction logic for MainWindow.xaml
28     /// </summary>
29     public partial class MainWindow : Window
30     {
31         AccelAmplifier amp;
32
33         //Vectors created for the required axis
34         List<double> xAccel = new List<double>();
35         List<double> yAccel = new List<double>();
36         List<double> zAccel = new List<double>();
37
38         List<double> addedAxis = new List<double>();//the result of the axis
39         added together
40
41         //Vectors created for the fake tremor experiment
42         List<double> noTremorX = new List<double>();
43         List<double> noTremorY = new List<double>();
44         List<double> noTremorZ = new List<double>();
45         List<double> lightTremorX = new List<double>();
46         List<double> lightTremorY = new List<double>();
47         List<double> lightTremorZ = new List<double>();
48         List<double> mediumTremorX = new List<double>();
49         List<double> mediumTremorY = new List<double>();
50         List<double> mediumTremorZ = new List<double>();
51         List<double> heavyTremorX = new List<double>();
52         List<double> heavyTremorY = new List<double>();
53         List<double> heavyTremorZ = new List<double>();
54
55         Timer t = new Timer(200);//timer created
```

```

55
56     /// <summary>
57     /// The main window for the program to run
58     /// </summary>
59     public MainWindow()
60     {
61         //Area for initializing variables
62         InitializeComponent();
63
64         /*
65         for (int i = 0; i < 2000; i++)// optional code to add zeros to the ↗
66             front of the vectors
67         {
68             xAccel.Add(0);
69             yAccel.Add(0);
70             zAccel.Add(0);
71         }*/
72
73         t.Elapsed += (o, e) =>
74         {
75             try {
76                 System.Windows.Application.Current.Dispatcher.Invoke(() =>
77                 {
78                     xPlotter.Plot(xAccel);
79                     yPlotter.Plot(yAccel);
80                     zPlotter.Plot(zAccel);
81                 });
82             } catch
83             {
84             }
85         };
86     }
87
88     /// <summary>
89     /// Actions to be carried out whenever the connect button is clicked. ↗
90     /// Searches for the accelerometer and amplifier and starts a timer.
91     /// </summary>
92     /// <param name="sender"></param>
93     /// <param name="eArgs"></param>
94     private void connectBTN_Click(object sender, RoutedEventArgs eArgs)
95     {
96         FrameDefinition.softwareEnvironment =
97             FrameDefinition.SoftwareEnvironments.CUSTOM;
98
99         amp = new AccelAmplifier();
100
101         amp.OnNewScopeData += Amp_OnNewScopeData;
102
103         amp.g_searchCompleted += (o, e) =>
104         {
105             var unfiltered = amp.getUnfilteredConnectedDevices();
106
107             if (unfiltered.Count == 1 && unfiltered[0] == "simulator")
108                 throw new Exception("no valid amplifier connected");
109         }
110     }

```

```

108         amp.SelectAmplifier(unfiltered[0]);
109
110         FrameDefinition.setForCSPWihtNautilus();
111         PatientManager.workingDirectory =           ↗
            Environment.CurrentDirectory;
112
113         amp.startAcquiring();//acquisition begins
114         t.Start();//timer start
115     };
116
117     amp.searchConnectedDevices();
118
119 }
120
121 private void Window_Closing(object sender,           ↗
    System.ComponentModel.CancelEventArgs e)
122 {
123     try
124     {
125         amp.stopAcquiring();
126     } catch{//closes the window whenever there is
127     {
128
129     }
130 }
131
132 /// <summary>
133 /// What happens when a new sample is read in from the amp
134 /// </summary>
135 /// <param name="sender"></param>
136 /// <param name="e"></param>
137 private void Amp_OnNewScopeData(object sender, newScopeDataArgs e)
138 {
139     // Optional code to remove the values at time=0 so that the graph ↗
140     // remains at the same timescale
141     //xAccel.Add(Math.Max(Math.Min(e.GetData()[0], 1), -1));
142     //xAccel.RemoveAt(0);
143     //yAccel.Add(Math.Max(Math.Min(e.GetData()[1], 1), -1));
144     //yAccel.RemoveAt(0);
145     //zAccel.Add(Math.Max(Math.Min(e.GetData()[2], 1), -1));
146     //zAccel.RemoveAt(0);
147
148     xAccel.Add(e.GetData()[0]);
149     yAccel.Add(e.GetData()[1]);
150     zAccel.Add(e.GetData()[2]);
151 }
152
153 /// <summary>
154 /// Actions to be carried out whenever the stop button is clicked.   ↗
155 /// Typically save to file.
156 /// </summary>
157 /// <param name="sender"></param>
158 /// <param name="e"></param>
159 private void stopBTN_Click(object sender, RoutedEventArgs e)
160 {

```

```

...c\acquisition project\accelerometer\MainWindow2.xaml.cs 4
160     string filename = Filenametb.Text+".csv";//input from text box.  ↗
        name for each experiment
161
162     string patientfolder = "patient" + Patientnumbertb.Text;// +  ↗
        @"\";//input from text box. patient number. used to create a new ↗
        directory
163
164     System.IO.Directory.CreateDirectory(@"C:\Users\Bens laptop  ↗
        \Documents\g.tec\acquisition project\" + patientfolder);// new  ↗
        directory created for each patient
165     amp.stopAcquiring();
166
167     //for (int i = 0; i < xAccel.Count; i++) // limit signal
168     //{
169         //    var s = xAccel[i];
170         //    if (s > 0.3) s = 10;
171         //    else s = -10;
172         //    xAccel[i] = s;
173     //}
174
175     for (int i = 0; i<50; i++) //Remove first 50 values to minimize  ↗
        the initial spike
176     {
177         xAccel.RemoveAt(i);
178         yAccel.RemoveAt(i);
179         zAccel.RemoveAt(i);
180     }
181
182     // Save to file within the specified directory for each patient.  ↗
        In the format date, time, axis, experiment name
183     Helper.Savefile(xAccel, @"C:\Users\Bens laptop\Documents\g.tec  ↗
        \acquisition project\" + patientfolder + @"\"+  ↗
        DateTime.Now.ToString("yyyyMMddHHmm") + "acc_x_" + filename);
184     Helper.Savefile(yAccel, @"C:\Users\Bens laptop\Documents\g.tec  ↗
        \acquisition project\" + patientfolder + @"\"+  ↗
        DateTime.Now.ToString("yyyyMMddHHmm") + "acc_y_" + filename);
185     Helper.Savefile(zAccel, @"C:\Users\Bens laptop\Documents\g.tec  ↗
        \acquisition project\" + patientfolder + @"\"+  ↗
        DateTime.Now.ToString("yyyyMMddHHmm") + "acc_z_" + filename);
186 }
187
188 /// <summary>
189 /// Offline post processing functions to be carried out whenever the  ↗
        Process button is clicked
190 /// </summary>
191 /// <param name="sender"></param>
192 /// <param name="e"></param>
193 private void process_Click(object sender, RoutedEventArgs e)
194 {
195     ///////////////////////////////////
196     //// Code for dealing with fake tremor experiment:  ↗
197
198     for (int patientnumber = 101; patientnumber < 106; patientnumber+  ↗
        +)//selection of patient number (101 to 105 represents  ↗
        experiments 1 to 5)

```

```

...c\acquisition project\accelerometer\MainWindow2.xaml.cs 5
199     {
200         string directory = @"C:\Users\Bens laptop\Documents\g.tec
\acquisition project\patient" + Convert.ToString
(patientnumber) + @"\";
201
202         List<double> scores = new List<double>();//create a list of
scores
203                                     //int i = 1000;
204
205         int i = 1000; // initial trimming point set to 1000
dataelements
206
207         //no tremor
208
209         string keyword = "notremor";
210
211         List<List<double>> chosenData = SearchAndLoad(directory,
keyword);//new list of lists containing the files returned
from the search and load function
212
213         List<double> notremor = Helper.AddAxis(Helper.Sum
(Helper.Abs(Helper.Window(chosenData[0])), true, i),
Helper.Sum(Helper.Abs(Helper.Window(chosenData[1])), true,
i), Helper.Sum(Helper.Abs(Helper.Window(chosenData[2])),
true, i));
214         //^^ chooses the X, Y and Z axis (indexed 0, 1, 2) of the
desired files based on a keyword. Processes each axis
moving ave, abs, sum. Adds each axis. Saves to csv.
Creates a new list.
215
216         double s_notremor = Helper.Score(notremor, false,
15000);//calculate scores from 15000 elements (60 secs at
fs = 250 Hz)
217
218         //light tremor
219
220         string keyword2 = "lighttremor";
221
222         List<List<double>> chosenData2 = SearchAndLoad(directory,
keyword2);//new list of lists containing the files
returned from the search and load function
223
224         List<double> lighttremor = Helper.AddAxis(Helper.Sum
(Helper.Abs(Helper.Window(chosenData2[0])), true, i),
Helper.Sum(Helper.Abs(Helper.Window(chosenData2[1])),
true, i), Helper.Sum(Helper.Abs(Helper.Window(chosenData2
[2])), true, i));
225         //^^ chooses the X, Y and Z axis (indexed 0, 1, 2) of the
desired files based on a keyword. Processes each axis
moving ave, abs, sum. Adds each axis. Saves to csv.
Creates a new list.
226
227         double s_lighttremor = Helper.Score(lighttremor, false,
15000);//calculate scores from 15000 elements (60 secs at
fs = 250 Hz)
228

```

```

229          //medium tremor
230
231          string keyword3 = "mediumtremor";
232
233          List<List<double>> chosenData3 = SearchAndLoad(directory, ➤
keyword3); //new list of lists containing the files ➤
returned from the search and load function
234
235          List<double> mediumtremor = Helper.AddAxis(Helper.Sum ➤
(Helper.Abs(Helper.Window(chosenData3[0])), true, i), ➤
Helper.Sum(Helper.Abs(Helper.Window(chosenData3[1])), ➤
true, i), Helper.Sum(Helper.Abs(Helper.Window(chosenData3 ➤
[2])), true, i));
236          //^^ chooses the X, Y and Z axis (indexed 0, 1, 2) of the ➤
desired files based on a keyword. Processes each axis ➤
moving ave, abs, sum. Adds each axis. Saves to csv. ➤
Creates a new list.
237
238          double s_mediumtremor = Helper.Score(mediumtremor, false, ➤
15000); //calculate scores from normally 15000 elements (30 ➤
secs at fs = 250 Hz)
239
240          //heavytremor
241
242          string keyword4 = "heavytremor";
243
244          List<List<double>> chosenData4 = SearchAndLoad(directory, ➤
keyword4); //new list of lists containing the files ➤
returned from the search and load function
245
246          List<double> heavytremor = Helper.AddAxis(Helper.Sum ➤
(Helper.Abs(Helper.Window(chosenData4[0])), true, i), ➤
Helper.Sum(Helper.Abs(Helper.Window(chosenData4[1])), ➤
true, i), Helper.Sum(Helper.Abs(Helper.Window(chosenData4 ➤
[2])), true, i));
247          //^^ chooses the X, Y and Z axis (indexed 0, 1, 2) of the ➤
desired files based on a keyword. Processes each axis ➤
moving ave, abs, sum. Adds each axis. Saves to csv. ➤
Creates a new list.
248
249          double s_heavytremor = Helper.Score(heavytremor, false, ➤
15000); //calculate scores from 15000 elements (60 secs at ➤
fs = 250 Hz)
250
251          //Scoring mechanism adds a score from each experiment to a ➤
new list of doubles "scores"
252
253          scores.Add(patientnumber);
254          scores.Add(i);
255          scores.Add(s_notremor);
256          scores.Add(s_lighttremor);
257          scores.Add(s_mediumtremor);
258          scores.Add(s_heavytremor);
259
260          Helper.Savefile(scores, directory + ➤
"trimmedscores1000.csv"); //save list of scores in a csv file ➤

```

```

        in the patient directory
261     }
262
263     ////////////////
264     //// Code for dealing with patient data
265
266     for (int patientnumber = 1; patientnumber < 6; patientnumber++)// 7
        selection of patient number
267     {
268
269         string directory = @"C:\Users\Bens laptop\Documents\g.tec 7
            \acquisition project\patient" + Convert.ToString 7
            (patientnumber) + @"\";
270
271         List<double> scores = new List<double>();//create a list of 7
            scores
272
273         // For loop which incrementally adjusts the initial trimming 7
            point of each signal thereby splitting the signals into 7
            subsections for internal analysis
274         for (int i = 1000; i <= 9000; i += 1000)//i = starting point 7
            for each measurement subsection
275         {
276             //Motion experiment - healthy arm
277
278             string keyword = "motionh";
279
280             List<List<double>> chosenData = SearchAndLoad(directory, 7
                keyword);//new list of lists containing the files returned 7
                from the search and load function
281
282             List<double> motionh = Helper.AddAxis(Helper.Sum 7
                (Helper.Abs(Helper.Window(chosenData[0])), true, i), 7
                Helper.Sum(Helper.Abs(Helper.Window(chosenData[1])), true, 7
                i), Helper.Sum(Helper.Abs(Helper.Window(chosenData[2])), 7
                true, i));
283             //^^ choses the X, Y and Z axis (indexed 0, 1, 2) of the 7
                desired files based on a keyword. Processes each axis 7
                moving ave, abs, sum. Adds each axis. Saves to csv. 7
                Creates a new list.
284
285             double s_motionh = Helper.Score(motionh, false, 1000);// 7
                calculate scores from x elements, typically set to 7500 7
                (30 secs at fs = 250 Hz)
286
287             //Motion experiment - affected arm
288
289             string keyword2 = "motiona";
290
291             List<List<double>> chosenData2 = SearchAndLoad(directory, 7
                keyword2);//new list of lists containing the files 7
                returned from the search and load function
292
293             List<double> motiona = Helper.AddAxis(Helper.Sum 7
                (Helper.Abs(Helper.Window(chosenData2[0])), true, i), 7
                Helper.Sum(Helper.Abs(Helper.Window(chosenData2[1])), 7

```

```

...c\acquisition project\accelerometer\MainWindow2.xaml.cs 8
true, i), Helper.Sum(Helper.Abs(Helper.Window(chosenData2
[2])), true, i));
294 //^^ chooses the X, Y and Z axis (indexed 0, 1, 2) of the
desired files based on a keyword. Processes each axis
moving ave, abs, sum. Adds each axis. Saves to csv.
Creates a new list.
295
296 double s_motiona = Helper.Score(motiona, false, 1000);//
calculate scores from x elements, typically set to 7500
(30 secs at fs = 250 Hz)
297
298 //Extension experiment - healthy arm
299
300 string keyword3 = "exth";
301
302 List<List<double>> chosenData3 = SearchAndLoad(directory,
keyword3);//new list of lists containing the files
returned from the search and load function
303
304 List<double> exth = Helper.AddAxis(Helper.Sum(Helper.Abs
(Helper.Window(chosenData3[0])), true, i), Helper.Sum
(Helper.Abs(Helper.Window(chosenData3[1])), true, i),
Helper.Sum(Helper.Abs(Helper.Window(chosenData3[2])),
true, i));
305 //^^ chooses the X, Y and Z axis (indexed 0, 1, 2) of the
desired files based on a keyword. Processes each axis
moving ave, abs, sum. Adds each axis. Saves to csv.
Creates a new list.
306
307 double s_exth = Helper.Score(exth, false, 1000);//
calculate scores from x elements, normally 2500 elements
(10 secs at fs = 250 Hz)
308
309 //Extension experiment - affected arm
310
311 string keyword4 = "exta";
312
313 List<List<double>> chosenData4 = SearchAndLoad(directory,
keyword4);//new list of lists containing the files
returned from the search and load function
314
315 List<double> exta = Helper.AddAxis(Helper.Sum(Helper.Abs
(Helper.Window(chosenData4[0])), true, i), Helper.Sum
(Helper.Abs(Helper.Window(chosenData4[1])), true, i),
Helper.Sum(Helper.Abs(Helper.Window(chosenData4[2])),
true, i));
316 //^^ chooses the X, Y and Z axis (indexed 0, 1, 2) of the
desired files based on a keyword. Processes each axis
moving ave, abs, sum. Adds each axis. Saves to csv.
Creates a new list.
317
318 double s_exta = Helper.Score(exta, false, 1000);//
calculate scores from x elements, normally 2500 elements
(10 secs at fs = 250 Hz)
319
320 //Scoring mechanism adds a score from each experiment to a

```



```

        new List<double>("scores")
321
322        scores.Add(patientnumber);
323        scores.Add(i);
324        scores.Add(s_motionh);
325        scores.Add(s_motiona);
326        scores.Add(s_exth);
327        scores.Add(s_exta);
328    }
329    Helper.Savefile(scores, directory + "trimmingscoresmot.csv");// save scores in a csv file in the patient directory
330
331    }
332
333    }
334    /// <summary>
335    /// Searches through a folder and finds csv files containing the keyword, loads them into vectors
336    /// </summary>
337    /// <param name="basepath"></param>filepath
338    /// <param name="keyword"></param>keyword for the search
339    /// <returns></returns>
340    private List<List<double>> SearchAndLoad(string basepath , string keyword)
341    {
342        var files = System.IO.Directory.GetFiles(basepath);//directory selection
343        var keywordFiles = files.Where(x => x.Contains(keyword)).ToList();//find files containing key words
344        keywordFiles.Sort();//sorted so that its in the order x, y, z - > 0, 1, 2
345        List<List<double>> keywordData = new List<List<double>>();//list of list of doubles
346        keywordData.Add(new List<double>());//one for each axis
347        keywordData.Add(new List<double>());
348        keywordData.Add(new List<double>());
349
350        //for (int fileIdx = 0; fileIdx < keywordFiles.Count(); fileIdx++)//.Count not working properly
351        for (int fileIdx = 0; fileIdx < 3; fileIdx++)//loop hard has been coded to 3 because there is always three axis
352        {
353            keywordData[fileIdx] = Helper.Loadfile(keywordFiles[fileIdx]);
354        }
355
356        return keywordData;
357    }
358
359    }
360    public class AccelAmplifier : Gtec2.AmpController//base station / amplifier communication functions. Sampling rate = 250
361    {
362        private volatile bool searchFinished = false;
363
364        public override GdsClientApiLibraryWrapper.GdsConfiguration createCustomConfiguration(ulong connectionHandle, string

```

```
        serialNumber, FrameDefinition.SoftwareEnvironments environment)
    {
365         GNautilusGdsClientApiLibraryWrapper.GdsGNutilusConfiguration
366         deviceConfiguration = new
367             GNautilusGdsClientApiLibraryWrapper.GdsGNutilusConfiguration();
368
369         deviceConfiguration.ValidationIndicator = false;
370         deviceConfiguration.BatteryLevel = false;
371         deviceConfiguration.AccelerationData = true;
372         deviceConfiguration.NumberOfScans = 0;
373         deviceConfiguration.SamplingRate = 250;
374
375         FrameDefinition.Fz = (int)deviceConfiguration.SamplingRate;
376         deviceConfiguration.Channels = new
377             GNautilusGdsClientApiLibraryWrapper.GdsGNutilusChannelConfigura
378             tion[GNautilusGdsClientApiLibraryWrapper.MaxNumberOfChannels];
379
380         for (int i = 0; i < deviceConfiguration.Channels.Length; i++)
381         {
382             deviceConfiguration.Channels[i].Enabled = false;
383             deviceConfiguration.Channels[i].Sensitivity =
384                 FrameDefinition.AmpSensitivity;
385             deviceConfiguration.Channels[i].BipolarChannel = -1;
386         }
387
388         //build and return complete GDS configuration structure
389         GdsClientApiLibraryWrapper.GdsConfiguration NAutilusConfig = new
390             GdsClientApiLibraryWrapper.GdsConfiguration();
391
392         NAutilusConfig.DeviceConfiguration = deviceConfiguration;
393         NAutilusConfig.DeviceInfo = new
394             GdsClientApiLibraryWrapper.DeviceInfo();
395         NAutilusConfig.DeviceInfo.DeviceType =
396             GdsClientApiLibraryWrapper.GdsDeviceType.GNautilus;
397         NAutilusConfig.DeviceInfo.Name = serialNumber;
398
399         return NAutilusConfig;
400     }
401 }
```